# Scale Two

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## ELECTRONIC PERSONAL SCALE

### OPERATING INSTRUCTIONS

Congratulations! You have just purchased fantastic Scale Two electronic personal scale.

Your new bathroom scale is extremely reliable and precise; it does not require any sort of maintenance.

### IMPORTANT SAFETY INSTRUCTIONS

Please read through the operating instructions carefully before you use the appliance for the first time and keep them safely for later reference; if necessary pass them on to the next owner.

- Handle your scale with care as it is a precision instrument. Therefore, avoid to expose it to extreme temperatures, damp, knocks, dust and direct sunlight. Do not drop the scale or knock on it.
- Keep the scale away from sources of heat.
- Do not jump on the scale.
- Take care that no liquids penetrate into the scale, as it is not waterproof. Do not leave or use the scale in very damp surroundings, humidity above 85% or contact with water and other liquid parts could cause corrosion and jeopardize the precision of the scale.
- Do not overload the scale.
- Do not open the scale. In the event of opening or tampering the guarantee becomes invalid.
- Go to the retailer for any maintenance.
- Keep the appliance out of the reach of children.
- Scale for domestic use.
- To improve the reading, weigh yourself twice, and if the weights are different, your weight is the middle of the two readings. Always weigh yourself on the same scale, at the same time of day, preferably undressed and before breakfast.
- Caution! The surface of the scale can be slippery, if wet.

### THE BODY WEIGHT

Checking the weight means performing an energetic balance: such balance results from the ratio between the quantity of ingested nutrients, i.e. the ingested calories, and the actually consumed calories. The daily necessary quantity of energy depends on the metabolism, on the life-style and on the type of activities performed by each single person. Controlling one's own weight is essential in order to protect health effectively: eating well means being healthy. It has been proved that most of diseases, many of which extremely dangerous, are associated to unbalanced eating habits not only in terms of quality but also in terms of quantity.

Among the most known, the following ones should be mentioned:

- ischemic heart diseases;
- cerebrovascular diseases;

- diabetes mellitus;
- intestine tumors;
- etc.

All these diseases are present above all among obese persons: therefore, controlling the weight together with controlling the quality of the nutrients becomes one of the essential factors to prevent diseases. Healthy eating means selecting the foodstuffs that supply the body with all the necessary nutriment to grow and keep the physiological functions and the biochemical processes balanced.

### HINTS FOR A HEALTHY FEEDING AND TO CONTROL WEIGHT

- To follow a balanced diet it is necessary not to exceed in the consumption of the same types of food neglecting the essential elements contained in other types of food, as well as not eating more than needed. A healthy diet includes suitable quantities of proteins, carbohydrates, fats, fibers, vitamins and minerals, as well as, obviously, water.
- Having a complete meal only once a day is a bad habit. It is better to have many small meals instead of just an abundant one. Three meals a day and a snack in the middle of the morning or of the afternoon represent the basis for a correct feeding.
- It is suggested to sit down and relax when eating. Select a comfortable environment can be relevant.
- Do not eat when you are focused on doing other activities, i.e. when you are watching TV or working, in order to avoid eating too much without even noticing it. To support the digestion process it is recommended to eat slowly and chew well, enjoying the color, taste and consistency of the food.
- Combining a correct feeding to physical activity is an essential condition to allow increasing the metabolism, consuming more calories and remaining in good shape.
- It is important to remember that with age human body seems to need less energy due to still not so clear reasons (maybe due to the fact that one leads a less active life): however, many aged people do not reduce their consumption of food and tend to grow fat.

Getting used to follow a balanced diet is surely of help in increasing the quality of life, promoting health and preventing diseases, as well as in increasing vitality and remaining in good shape!

### CAUTION!

- Never follow diets autonomously, always refer to a doctor or dietician.
- Self-measurement means control, not diagnosis or treatment. Unusual values have always to be discussed with one's own physician. Under no circumstances must you change the dose of medicines prescribed by your Doctor.

## OPERATION

1. Select the weighing unit of measure (kg, lb or st) with the key located on the bottom of the scale. It is possible to select the unit of measure if the scale is turned off as well as it is turned on.

2. Place the scale on a flat, hard surface. If the floor is loose or uneven (e.g. mats, carpets, linoleum) the exactness of your weight is not guaranteed.

3. Press the centre of the scale to turn it on. On the display screen appears "CAL".

4. Wait a few seconds and when "0" (or "0.0") appears on the display screen, the scale is ready for use.

5. Step onto the scale with your feet parallel and your weight equally distributed (see illustrations), before the zeros disappear. Stand still while the scale measures your weight.

Do not lean on anything while being weighed as this could alter the reading.

6. Step off the scale; the display screen will continue to show your weight for a few seconds then the scale will turn itself off automatically.

### WARNING

- When the display screen shows "LO" or the symbol of the battery, it's necessary to replace the batteries bearing in mind the polarity indicated (+) / (-).
- In case of overweight, the display screen shows "O-Ld", remove the weight from the scale to avoid any damage of the mechanism.
- When the display screen shows "Err", step off the scale, wait a few seconds and weight again.

## BATTERY INSTALLATION / REPLACEMENT

This scale operates with replaceable batteries.

1. Insert three 1,5 V AAA alkaline batteries in the battery compartment at the bottom of the scale, bearing in mind the polarity indicated.

2. When the batteries are discharged the display will show either "LO" or the battery symbol. Remove the batteries if the scale are not to be used for a prolonged period of time. Always dispose of the batteries in the special containers.

## CLEANING

- Clean the platform with a damp cloth.
- Clean the bottom of the scale with just a soft dry cloth.
- To clean the scale never use chemical or abrasive products.

## SPECIFICATIONS

3 x 1,5 volt AAA alkaline batteries

Reverse LCD display

Maximum tolerance: +/-1% + 1 division

Working conditions: +10°C +40°C; RH 85%

Storage conditions: +10°C +50°C; RH 85%

All rights for any technical modifications are reserved

The conformity, certified by the CE-marking on the device, refers to the EEC Directive 2004/108 on electromagnetic compatibility.

## REPAIRS/DISPOSAL

- Repairs to electrical appliances must only be carried out by a qualified electrical technician. If improper repairs are carried out, the guarantee lapses and all liability is refused.
- Never operate the appliance after it has malfunctioned, if it has been dropped or has been damaged in any other way (cracks/breaks in the casing).
- Do not push any object into the appliance, do not disassemble the appliance.
- If the appliance is beyond repair, make it unusable immediately and hand it in at the collection point provided for the purpose.
- Do not dispose of electrical appliances in the domestic waste. Please use the communal collection points.
- Ask your local authority if you do not know where the communal collection points are located.
- If electrical appliances are disposed of thoughtlessly, exposure to the elements may cause hazardous substances to find their way into the ground-water and from there into the food chain, or to poison the flora and fauna for many years to come.
- If you are replacing the appliance with a new one, the seller is legally obliged to take back the old one at least free of charge for disposal.

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Making life colorful and enjoy is the goal and endeavour of us all. This also includes designer items of practical use: This is our guiding principle!

A lot of household appliances are important and useful occupants of our homes however, they are often the ugly ducklings of the family. That's why we came up with the idea to develop attractive and appealing household products.

All our products are invented and engineered in Switzerland

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